



Monofocal IOLs

Aspheric HP IOLs

* SeeLens HP



* BunnyLens HP



* Available in preloaded configuration

	SeeLens HP *	BunnyLens HP *
Platform	C-Loop	4 - Loop
Overall Diameter	13 mm	11 mm (>10D) 11.5mm (≤10D)
Placement	Capsular Bag	
Optic Diameter	6 mm	
Power Range	+5 to +10 (1 increments) +10 to +30 (0.5 increments) +30 to +35 (1 increments)	
Optic Design	Aspheric	
360° Continuous Square Edge	Yes	
Haptic Angulation	5°	
Material	Hydrophobic acrylic with bonded UV absorber and violet light filter	
Refractive Index	1.48 (@ 35°C)	
A-constant (SRK/T) for Optical or Immersion US Biometry	119	118.9
A-constant (SRK/T) for Contact US Biometry	118.5	118.4
Sterilization	ETO	


Providing Excellent Vision Quality For Your Patient

The Aspheric HP IOLs were designed by the most advanced optical tools to provide:

- reduced spherical aberration
- functional vision in night conditions

The Aspheric Hydrophobic lens from Hanita Lenses provides the patient with an excellent vision quality at day and night conditions by using state-of-the-art aberration-free aspherical optical design.

Stability and Centration



PRELOADED HYDROPHOBIC IOL SYSTEM

- Simplifying cataract surgery
- Ready for implant
- Smooth and EASY injection through a 2.2 mm incision