# Accessories

### AssiAnchor

#### **Attributes**

- Hanita Lenses upgraded capsular anchor device
- Securing the capsular bag to the scleral wall
- Enables cataract surgery for subluxated lens
- Centers subluxated capsule
- Provides a wide contact between the device and the anterior capsule
- · CE approved

# Capsule, Sulcus<sup>1,2</sup> and AssiAnchor dimensions

Capsulorhexis


Capsule

Ciliary body

### **Specifications**

Material: Blue tinted PMMA YAG laser: compatible

Developed with Prof. E.Assia, Israel



# **Upgraded Capsular Anchor Device**







#### Ref 1,2

- Assia E.I., Apple D.J.: Side-View Analysis of the Lens.
  Part 1: The Crystalline Lens and the Evacuated Bag Arch. Ophthalmol. 1992; 110: 89-93.
- Assia E.I., Apple D.J.: Side-View Analysis of the Lens. Part 2: Positioning of Intraocular Lenses. Arch. Ophthalmol. 1992; 1 10: 94-98.

